

The Ramp

ORIF of Tibia, Ankle, Midfoot and Forefoot FX

Save time with The Ramp while positioning and imaging the operative leg for a variety of fractures.

- Consistent and fast positioning every time
- Used for ORIF of tibia, ankle, mid-foot and forefoot fractures
- Variety of options including Adult, Pediatric, Prone Positioner, and Surgery Center ramp
- Can be used with Bone Foam Custom Sterile Bags

“The leg elevator helps me save time by positioning and stabilizing the operative leg above the nonoperative leg so lateral C-Arm Imaging is consistent and fast.”



The Pediatric Ramp

Orif Of Tibia, Ankle, Midfoot & Forefoot FX

Pediatric Ramp
#925.PTR

Pediatric Surgery
Center Ramp
#925.PSCTR



Adult Ramp #925.TR

Surgery Center Ramp #925.SCTR

Ramp Shorty #925.TRS



Custom Sterile Bags

Now Available



The Ramp Prone Positioner

Orif Of Tibia, Distal Tibia, Ankle

The Ramp
Prone Positioner
#925.TRPP

The difference between the Adult and Surgery Center sizes, is the Surgery Center is designed for facilities that use the mini c-arm, and is 2 inches lower in height compared to the Adult size. The Adult size is used in most hospitals.