


NEW
PRODUCT



Mayo Band™

The Mayo Band™ reduces the risk of pressure on the soft tissue and helps avoid skin and nerve damage when the extremities are positioned on the hard Mayo Stand platform during orthopedic procedures.

- Soft tissue friendly, made of memory foam.
- Easy to use! Simply place the pad on the stand.
- Sterile Bag available for a simple set up

Custom Sterile Bags



Mayo Band™ Sterile Bags
#925.MBSB



Designed
with a
SPEEDBUMP
to hold instruments
during procedure

PATENT PENDING
Mayo Stand not included

Transform the Mayo Stand into a safe, soft-tissue friendly surface for positioning extremities.

- Applicable across numerous orthopedic surgeries: hip or arm abduction in ORIF acetabulum or total hips, or for ORIF of upper extremity conditions
- Various sizes to ensure proper fit with your existing Mayo Stand

Mayo Band™ Large #925.MBL

Mayo Band Large fits mayo stand trays that measure 20"x25"

Mayo Band™ Regular #925.MBR

Mayo Band Regular fits mayo stand trays that measure 16"x21" up to 16.75"x21.5"

Mayo Band™ Small #925.MBS

Mayo Band Small fits mayo stand trays that measure 12"x19" up to 13"x19.625"



Mayo Band Sterile Bags
(Case of 20, fits all
Mayo Band sizes)
#925.MBSB

Mayo Stand not included